



Technical Note DOSE 2_TN001_151124 01

Subject: DOSE² RS232 Command Set
Product affected DOSE²
Product functions affected Communication with a PC via RS232 port
Type of action: Follow the instruction in this Technical Note
Technical Note release: 2015-11-24, Schwarzenbruck

The DOSE² electrometer features the ability to exchange commands and communication with a PC via an RS232 port. When using applications such as Tera Term Pro to facilitate communication between the host PC and DOSE², the PC RS232 interface should be set to 19200 baud, 8 databits and no parity. All commands are ASCII based.

Command Structure

The command structure is defined as a stream of bytes with the structure:

Start Char	Command	Stop Char
------------	---------	-----------

Where:

Start Char '<'

Command see *Command Summary* section

Stop Char '>'

Response Structure

All responses have the same structure called the Response Structure. The Response Structure is defined as:

Command Status Character	Command Response
--------------------------	------------------

Where:

Status Character Indicates whether previous command was successful.

Status Character	Description
*	Command received and executed normally.
!	Command received but could not execute
?	Unknown Command

Command Response Command result if applicable, see *Command Detail* section.



Command Summary

1)	<GID>	Get the identification of the electrometer
2)	<GSN>	Get the serial number of the electrometer
3)	<GRGx>	Get the selected range for a specified channel
4)	<GBSx>	Get the selected bias setting for a specified channel
5)	<GBVx>	Get the selected measured bias value for a specified channel
6)	<GRx>	Get the rate for a specified channel
7)	<GCx>	Get the charge for a specified channel
8)	<GDRx>	Get the doserate for a specified channel
9)	<GDx>	Get the dose for a specified channel
10)	<GCS>	Get the charge status
11)	<GCT>	Get charge type
12)	<GZS>	Get the zeroing status
13)	<SRGx y>	Set input range for specified channel
14)	<SBSx y>	Set bias setting for specified channel
15)	<SCT x>	Set charge collection type
16)	<SVx>	Set channel view
17)	<DZ>	Do zeroing
18)	<STRC>	Start charge collection
19)	<STPC>	Stop charge collection
20)	<EDCx>	Data recording toggle
21)	<SRUx>	Unfiltered rate output

Command Detail

1. Get Identification

Get the identification of the electrometer

Syntax: <GID>

Parameters: None

Return: "DOSE2"

Example: <GID>*DOSE2

2. Get Serial Number

Get the serial number of the electrometer

Syntax: <GSN>

Parameters: None

Return: 7 digit serial number

Example: <GSN>*0123456

3. Get Range

Get the selected range for a specified channel

Syntax: <GRGx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: H high range

L low range

Example: <GRG1>*H



4. Get Bias Setting

Get the selected bias setting for a specified channel

Syntax: <GBSx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: bias setting in volts

Example: <GBS1>*150

5. Get Bias Value

Get the selected measured bias value for a specified channel

Syntax: <GBVx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: bias value in volts

Example: <GBV1>*-152

6. Get Rate

Get the rate for a specified channel

Syntax: <GRx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: measured rate

Example: <GR1>*-0.011 nA

7. Get Charge

Get the charge for a specified channel

Syntax: <GCx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: measured charge

Example: <GC1>*-0.082 nC

8. Get DoseRate

Get the dose rate for a specified channel

.

Syntax: <GDRx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: measured doserate

Example: <GDR1>*-0.082 Rm²/hA

9. Get Dose

Get the dose for a specified channel.

Syntax: <GDx>

Parameters: x: 1 = Channel 1, 2 = Channel 2

Return: measured dose

Example: <GD1>*-0.082 Ci



10. Get Charge Status

Get the charge status

Syntax: <GCS>

Parameters: none

Return: charge status

Status Character	Description
I	Idle
C	Collection in progress
A	Armed and waiting for trigger

Example: <GCS>*I

11. Get Charge Type

Get charge type

Syntax: <GCT>

Parameters: none

Return: charge type

T x timed, x = duration in seconds

C Continuous

TRG trigger

Example: <GCT>*T 15

12. Get Zero Status

Get the zeroing status

Syntax: <GZS>

Parameters: none

Return: zeroing status

0 = zeroing not in progress

1 = zeroing in progress

Example: <GZS>*1

13. Set Range

Set input range for specified channel

Syntax: <SRGx y>

Parameters: x: channel, 1 = channel 1, 2 = channel 2

y: range, H = high range, L = low range

Return: status (see response structure)

Example: <SRG1 L>*

14. Set Bias Setting

Set bias setting for specified channel

Syntax: <SBSx y>

Parameters: x: channel, 1 = channel 1, 2 = channel 2

y: bias setting in volts, +-1000

Return: status (see response structure)

Example: <SBS1 150>*



15. Set Charge Type

Set charge collection type

Syntax: <SCT x>

Parameters: x: charge type C = continuous

Ty = timed, y = duration in seconds

TRG = trigger

Return: status (see response structure)

Example: <SCT T15>*

16. Set Channel View

Set channel view

Syntax: <SVx>

Parameters: x: channel view

1 = Channel 1view

2 = Channel 2 view

B = dual channel view

Return: status (see response structure)

Example: <SVB>*

17. Do Zero

Do zeroing

Syntax: <DZ>

Parameters: none

Return: status (see response structure)

Example: <DZ>*

18. Start Charge Collection

Start charge collection

Syntax: <STRC>

Parameters: none

Return: status (see response structure)

Example: <STRC>*

*Note: with DOSE2 software version 2.0, a function was added which requires pressing of the **Clear** button after the completion of a charge collection before starting a new collection. From the command side, this is accomplished by sending an additional <STRC> command.*

19. Stop Charge Collection

Stop charge collection

Syntax: <STPC>

Parameters: none

Return: status (see response structure)

Example: <STPC>*



20. Data Recording Toggle

Enable or disable recording data to the Add List when using the standard measurement Mode

Syntax: <EDCx>

Parameters: x: enable/disable

0 = Disable data recording

1 = Enable data recording

Return: status (see response structure)

Example: <EDC0>*

Note: When the DOSE2 is turned off and back on via the power button or via unplugging re-plugging, the data recording is automatically re-enabled.

21. Unfiltered Rate Output

Enable or disable streaming display of unfiltered (10 Hz) rate data for both channels

Syntax: <SRUx>

Parameters: x: enable/disable

0 = Disable streaming

1 = Enable streaming

Return: channel1value,channel2value (values shown in fA)

Example: <SRU1>852,-1653

In case of any further questions, please contact us under:

IBA Dosimetry GmbH
Service Department
Bahnhofstrasse 5
DE-90592 Schwarzenbruck, Germany
Phone: +49 9128 607 38
Fax: + 49 9128 607 26
Mailto: service@iba-group.com