



Technical Note DOSE1_TN002_200428 01

Subject: DOSE1 Software compatibility with Windows 10 (64 bit)
Product affected: DOSE1
Function affected: -
Type of action: Follow the instruction in this Technote

For the computer system requirements of the DOSE 1 software, it has been tested on Windows 10, 64 bit operating system (OS). The compatibility to this OS has been approved if following requirements are fulfilled:

1. Ensure that .NET Framework 1.1 is installed and running.
2. Start the DOSE1 software as described in the User's Guide.

For any further questions, please contact us:

IBA Dosimetry GmbH
Service Department
Bahnhofstrasse 5
90592 Schwarzenbruck, Germany
Tel.: +49 9128 607 38
Email: service-emea@iba-group.com